



Disinfecting Procedure for

ACCUSCAN DP™

Effective Date: 10/25/23

Purpose of Policy:

To ensure that appropriate cleaning is performed to the Accuscan DP iPad and Scanning Device to aid in preventing the transmission of infection, and to ensure the equipment is maintained in optimal functionality and cleanliness.


Policy:

Disinfecting Products & Equipment Needed:

1. 70% isopropyl alcohol wipes

Procedure:

1. After Patient Use Perform Hand Hygiene
2. Disinfect the Accuscan DP iPad and Scanning Device using 70% isopropyl alcohol wipe, gently wipe the hard, non-porous surfaces
3. Place disinfected Accuscan DP iPad and Scanning Device in secure location for future use



Sergio Ramos

Co-Owner & CTO

3D Orthotics